

Oven

Constant-Temperature Drying Oven

Features:

- * External material: cold-rolled steel with anti-bacteria powder coating.
- * Stainless steel inner chamber, round angle structure, adjustable shelves.
- * PID control with LED display, soft touch operational key, timing function and over-temperature protection.
- * Door with double-layer tempered glass observation window, adjustable air-tight buckle lock, and silicon rubber seal, which provide good seal ability and efficiently prevent heat loss.



Technical Parameters:

Model	BOV-T30C	BOV-T70C	BOV-T140C	BOV-T270C
Capacity	43L	70L	137L	270L
Temp. Range	50~200°C			
Temp. Fluctuation	±1°C			
Ambient Temp.	5~40°C			
Timing Range	1~9999min			
Shelves No.	2 pcs		3 pcs	
Consumption	850W	1000W	2000W	3000W
Power Supply	Standard: 220V±10%, 50/60Hz; Optional: 110V±10%, 60Hz			
Internal Size(W*D*H) mm	350*350*350	400*390*450	500*500*550	600*600*750
External Size(W*D*H) mm	635*615*552	680*653*631	830*665*772	930*765*972
Package Size(W*D*H) mm	780*750*700	870*700*800	960*780*950	1060*880*1150
Gross Weight(kg)	60	70	105	145

Forced Air Drying Oven

Features:

- * External material: cold-rolled steel with anti-bacteria powder coating.
- * Stainless steel inner chamber, round angle structure, adjustable shelves.
- * Double-layer glass observation window, adjustable airtight buckle lock.
- * PID control with LED display.
- * Over-temperature protection.



Technical Parameters:

Model	BOV-T25F	BOV-T50F	BOV-T105F	BOV-T200F
Capacity	25L	54L	108L	212L
Temp. Range	50~200°C			
Temp. Precision	0.1°C			
Temp. Fluctuation	±1°C			
Ambient Temp.	5~40°C			
Timing Range	1~9999min			
Shelves No.	2 pcs			
Consumption	650W	1000W	1600W	2500W
Power Supply	Standard: 220V±10%,50/60Hz; Optional: 110V,60Hz			
Internal Size(W*D*H) mm	298*290*274	419*332*348	550*348*550	602*498*602
External Size(W*D*H) mm	580*554*455	708*595*532	830*620*735	880*775*785
Package Size(W*D*H) mm	710*660*620	840*690*700	960*720*900	1010*880*940
Gross Weight(kg)	51	66	115	145