

Pharmacy Refrigerators HYC-1378




Lockable Design



1378L Large Space Storage



Superior Cabinet Temperature
Uniformity, Auto Defrost



The test point of performance in the cabinet
increased from 2 points to 10 points



LED Temperature Display,
Display Precision 0.1°C



Multiple Alarm Functions :

Including high temperature alarm, low temperature alarm, remote
remote alarm, sensor error alarm, power failure alarm, low battery
power alarm, door ajar alarm.



Energy Efficient Hydrocarbon (HC)
Refrigeration System

Optimised refrigeration system, utilising high-efficiency HC
compressor significantly reduces energy consumption.



Greener by Design

HC refrigerants are natural non-toxic and more environmentally
friendly than HFC refrigerants.

Specifications

HYC-1378

Technical Data

Cabinet Type	Upright, double door	
Climate Class	N	
Cooling Type	Forced air cooling	
Defrost Mode	Auto	
Refrigerant	HC	
	220v	115v
Noise (dB (A))	43	45

Performance

Temp Range (°C)	2~8	
-----------------	-----	--

Control

Controller	Microprocessor	
Display	LED	

Electrical Data

Power Supply (V/Hz)	220~240V/50&60Hz	115/60
Power (W)	600	850
Electrical Current (A)	3	7.35

Dimensions

Capacity (L/Cu.Ft)	1378/48.7	
Net/Gross Weight (approx)	310/360 (kg) 683.4/793.7 (lbs)	
Interior Dimensions (W*D*H)	1320*700*1500 (mm) 52.0*27.6*59.1 (in)	
Exterior Dimensions (W*D*H)	1440*950*1980 (mm) 56.7*37.4*78.0 (in)	
Packing Dimensions (W*D*H)	1570*1010*2120 (mm) 61.8*39.8*83.5 (in)	
Container Load (20'/40'/40'H)	7/14/14	



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Porthole	Y
Shelves /Basket	12/-
Temperature Recorder	Optional
Printer	Optional
USB Interface	Y

Other

Certificate	CE	UL
-------------	----	----

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International